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Logowanie 1-2-3

The image shows the LabScala login interface. At the top is a blue header with the LabScala logo and name. Below the header, the text "Login with username and password" is displayed. There are two input fields: "Username:" and "Password:". Below the "Password:" field is a blue "Login" button and a link "Forgot your password?".

AUTH.01.1

The image shows the LabScala password change interface. At the top is a blue header with the LabScala logo and name. Below the header, the text "Password change" is displayed. There are three input fields: "Old password:", "New password:", and "Repeat new password:". Below the "Repeat new password:" field are two buttons: "Change" and "Back".

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My Schemes My registry Administration

Welcome to the new LabScala user interface

My EQA

Links & notes

LATEST 3 NOTIFICATIONS

You have no notifications.

NEXT 3 CLOSING SCHEMES

Round	Closing date	Status
Clinical chemistry: Endocrinology Hormones A, June, 2012	30.12.2012	Draft
Clinical chemistry: Endocrinology Hormones B, June, 2012	30.12.2012	Draft

MY REPORTS

Place orders Fill results

View reports Go to Mainio

IMPORTANT NOTES

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
Additional info

Labquality thanks all the congress participants and exhibition guests for visiting our booth and presentations at BALM 2012!

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My Schemes My registry Administration

Devices list

Help

▼ Devices

New device

Display: Active ▼

Nickname	Instrument	Manufacturer	Products	Client	
<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>

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[-] Hormones A

Display: All

Analyte	Reagent	Calibrator	Measurement unit	Active
Digoxin	Abbott Architect		nmol/l	<input checked="" type="checkbox"/>
Ferritin				<input checked="" type="checkbox"/>
Folate				<input checked="" type="checkbox"/>
hCG total				<input checked="" type="checkbox"/>
hCG intact				<input checked="" type="checkbox"/>
T3				<input checked="" type="checkbox"/>
T3 free				<input checked="" type="checkbox"/>
Thyroxine	Abbott Architect		nmol/l	<input checked="" type="checkbox"/>
Thyroxine, free	Abbott Architect		pmol/l	<input checked="" type="checkbox"/>
TSH	Abbott Architect		mU/l	<input checked="" type="checkbox"/>

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My Schemes My registry Administration

Add device's reagent info

Help

▼ Add device's reagent info for Digoxin

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DEVICE

Manufacturer: Architect i2000SR
Instrument: ABBOTT
* Nickname: Archie

* Start date: 17.07.2012
End date:
Location: Room 345
Client system ID: 1234567

USED FOR MEASURING SCHEMES

Product	Active
[+] Hormones A	<input checked="" type="checkbox"/>
Hormones B	<input type="checkbox"/>

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Analizator wybrany

The screenshot shows a web application interface with a navigation bar at the top containing 'My Schemes', 'My registry', and 'Administration'. Below the navigation bar is the 'Devices list' section, which includes a 'Help' link. A blue header bar for the 'Devices -' section contains a 'New device' button. Below this, a status bar indicates 'Showing entries 1-1 [Total 1]' and a 'Display: Active' dropdown menu. The main content is a table with the following columns: Nickname, Instrument, Manufacturer, Products, and Client. A search icon is located to the right of the table. Below the table, another 'New device' button is visible.

Nickname	Instrument	Manufacturer	Products	Client
Archie	Architect i2000SR	ABBOTT	Hormones A	

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
The screenshot shows the LabScala user interface. At the top, there is a navigation bar with a home icon, 'My Schemes', 'My registry', and 'Administration'. Below this, a dropdown menu for 'My Schemes' is open, showing options: 'Fill results', 'View all', and 'Mainio'. The main content area is titled 'LabScala user interface' and features a section for 'My EQA'. Under 'My EQA', there are two sub-sections: 'LATEST 3 NOTIFICATIONS' (with the message 'You have no notifications.') and 'NEXT 3 CLOSING SCHEMES'. The 'NEXT 3 CLOSING SCHEMES' section contains a table with the following data:

Round	Closing date	Status
Clinical chemistry: Endocrinology Hormones A, June, 2012	30.12.2012	Draft

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Welcome to the Clinical Chemistry Endocri

Request>>Instructions>>Pre-analytics>>Analytics>>Post-analytics>>Exit

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Welcome to the Clinical Chemistry Endocrinology Survey

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MY SCHEME SPECIFIC CONTACTS

Name	Send E-mail notification to	
Result filler	results@lab.com	 

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Clinical chemistry: Endocrinology Instru

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GENERAL INFO

Welcome to the hormones A scheme of August! If the kit is incomplete,

SAFETY INFO

Quality control specimens derived from human blood must be handled
are found to be HBs-Ag and HIV-Ab negative when tested with license
these or other infectious diseases.

SCHEDULE

Analyze the samples as soon as possible and report the results no late

SAMPLE INFORMATION

Samples A1 and A2 are lyophilised human sera. The unopened lyophilis
and A2. Replace the stopper and allow the bottles to stand for 30 minut
batches. Freeze immediately. A little before analysing, thaw, mix and al
stable for at least two weeks at -20 °C.

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Clinical chemistry: Endocrinology Preanalytics

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Sample registration

▼ Sample registration

Product	Code	Quantity ordered	Sample arrival date	Quantity received	Sample storage conditions	Sample preparation date
Hormones A, June 2012	2300	5	16.07.2012	7 <input type="text"/>	Refrigerator	16.07.2012

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
First Previous 1 Next Last

Sample_A1 Sample_A2

▼ Hormones A measurement results for Sample_A1

Measurement date	Analyte	Device	Instrument	Reagent	Result	Unit
	Digoxin	--Choose--				
		--Choose--				
		--Choose--				
		--Choose--				
		--Choose--				
		--Choose--				
		--Choose--				
		--Choose--				

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The screenshot shows a web form titled "Add a device". At the top, there is a blue header with a dropdown arrow and the text "Add a device". Below this, the "Client" section includes a "Client code:" label and an input field with search and refresh icons. The "DEVICE" section contains several fields: "Instrument:" and "Manufacturer:" (both with dropdown icons), "Nickname:", "Start date:" (with a calendar icon), "End date:" (with a calendar icon), "Location:", and "Client system ID:". Below the device fields is a section titled "USED FOR MEASURING SCHEMES" which contains a "Product" dropdown menu currently set to "Hormones A" and an "Active" checkbox with a plus icon.

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USED FOR MEASURING SCHEMES

Product	Active
[+] Hormones A	<input checked="" type="checkbox"/>
Hormones B <input type="text"/>	<input type="checkbox"/>

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The screenshot shows a web application interface with a blue header bar containing three tabs: "My Schemes", "My registry", and "Administration". Below the header, the main content area is titled "Add device`s reagent info" with a "Help" icon on the right. A blue bar below the title indicates the current test: "Add device`s reagent info for Digoxin". The form contains the following fields:

- * Reagent manufacturer: A text input field containing "ABBOTT" with a clear button (X) and a dropdown arrow.
- * Reagent: A text input field containing "Abbott Architect" with a clear button (X) and a dropdown arrow.
- * Measure unit: A dropdown menu currently showing "nmol/l".
- Methodics: A text input field containing "Abbott Architect".

At the bottom of the form, there are two buttons: "Back" on the left and "Save" on the right.

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Thyroxine	Abbott Architect		nmol/l	<input checked="" type="checkbox"/>
Thyroxine, free	Abbott Architect		pmol/l	<input checked="" type="checkbox"/>
TSH	Abbott Architect		mU/l	<input checked="" type="checkbox"/>

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Clinical chemistry: Endocrinology: Hormones B June, 2012 Validation

[Help](#)

Request>>Instructions>>Pre-analytics>>Analytics>>Post-analytics>>Exit

Validate results

Validation results

Product	Sample set	Sample	Errors	Last saved	Date	Status
Preanalytics		Registration	OK	Superadmin, Client	17.07.2012 07:04	Accepted
Hormones B	1	Sample_B1	OK	Superadmin, Client	17.07.2012 16:47	Draft
Hormones B	1	Sample_B2	Row "Aldosterone" result is out of limit boundaries. Please check result and unit.	Superadmin, Client	17.07.2012 16:47	Error

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Validate results

Validation results

Product	Sample set	Sample	Errors	Last saved	Date	Status
Preanalytics		Registration	OK	Superadmin, Client	17.07.2012 07:04	Accepted
Hormones B	1	Sample_B1	OK	Superadmin, Client	17.07.2012 12:13	Accepted
Hormones B	1	Sample_B2	OK	Superadmin, Client	17.07.2012 12:14	Accepted

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The screenshot displays the LabScala user interface. At the top, there is a navigation bar with a home icon and three tabs: "My Schemes", "My registry", and "Administration". Below the navigation bar, a welcome message reads "Welcome to the new LabScala user interface" with a "Help" link on the right. The main content area is divided into two panels. The left panel, titled "My EQA", contains sections for "LATEST 3 NOTIFICATIONS" (stating "You have no notifications."), "NEXT 3 CLOSING SCHEMES" (a table with 3 rows), and "MY REPORTS" (a link to "Go to Mainio reports"). The right panel, titled "Links & notes", contains sections for "SHORTCUTS" (links for "Place orders", "View reports", "Fill results", and "Go to Mainio"), "IMPORTANT NOTES" (three notification items), and a footer link for "Additional info".

My Schemes My registry Administration

Welcome to the new LabScala user interface Help

My EQA

LATEST 3 NOTIFICATIONS

You have no notifications.

NEXT 3 CLOSING SCHEMES

Round	Closing date	Status
Clinical chemistry: Endocrinology Hormones A, June, 2012	30.12.2012	Draft
Clinical chemistry: Endocrinology Hormones B, June, 2012	30.12.2012	Draft

MY REPORTS

Go to Mainio reports

Links & notes

SHORTCUTS


- Place orders
- View reports
- Fill results
- Go to Mainio

IMPORTANT NOTES

- NEW IT-SYSTEM: LabScala**
Dear client,
Labquality will took into use its new IT-system, LabScala, in May 2012. In the first phase of LabScala you will have access to some electronic forms as well as user administration. You will also have a manual for LabScala use accessible from the top right corner of every page, under the "Help" – text. Please do not hesitate to contact us if there are any questions concerning the use of LabScala. We hope to be able to provide you with an even better service with the help of LabScala.
- Scheme 5260, Fungal culture 2, 2012, delivery postponed to June.**
Fungal culture 2, 2012 (5260, May) delivery has been postponed to June. The new delivery date will be June 12, 2012.
- BALM2012 succesful**
Labquality thanks all the congress participants and exhibition guests for visiting our booth and presentations at BALM 2012!

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